



## User's Manual

# VIM IgG (Human) ELISA Kit



DEIA4437



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

Anti-VIM IgG (Human) ELISA Kit is an ELISA kit for a semiquantitative detection of IgG antibodies to the intermediate filament protein vimentin in human serum.

### Storage

Store the kit and the kit reagents at 2-10°C, in a dry place and protected from the light. For more detailed information, please download the following document on our website.

### Precision

Intra-assay Precision: 2.8%-4.4%

Inter-assay Precision: 4.3%-11.3%