



## User's Manual

# HIA IgG ELISA Kit



DEIA461



96T



This product is for research use only and is not intended for diagnostic use.



For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

HIA IgG ELISA Kit, is a qualitative assay intended for the detection of heparin dependent antibodies of the IgG isotype, in human plasma, by clinical laboratories. It is intended for in vitro diagnostic use.

### Storage

Store the kit at 2-8°C. For more detailed information, please download the following document on our website.