



User's Manual

Mouse Acenaphthene IgG ELISA Kit

REF

DEIA478



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

Mouse Acenaphthene IgG ELISA Kit is an ELISA kit intended for a qualitative detection of anti-Acenaphthene IgG antibodies in mouse serum. The kit is for research use only.

Storage

Store the kit reagents at 2-10°C, in a dry place and protected from the light. Avoid freezing.