



User's Manual

Mouse Phenanthrene IgG ELISA Kit

REF

DEIA480



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

Mouse Phenanthrene IgG ELISA Kit is an ELISA kit intended for a qualitative detection of anti-Phenanthrene IgG antibodies in mouse serum.

Storage

Store the kit reagents at 2-10°C, in a dry place and protected from the light. Avoid freezing.