



User's Manual

Mouse Fluoranthene IgG ELISA Kit



DEIA481



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

Mouse Fluoranthene IgG ELISA Kit is an ELISA kit intended for a qualitative detection of anti-Fluoranthene IgG antibodies in mouse serum.

Storage

Store the kit reagents at 2-10°C, in a dry place and protected from the light. Avoid freezing.