



User's Manual

Thromboxane B2 Express ELISA Kit



DEIA5035



96T, 480T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

This assay is based on the competition between TXB2 and a TXB2-acetylcholinesterase (AChE) conjugate (TXB2 Tracer) for a limited number of TXB2 Monoclonal Antibody.

Storage

This kit will perform as specified if stored as directed at -20°C and used before the expiration date indicated on the outside of the box. For more detailed information, please download the following document on our website.

Detection Range

15.6-2000 pg/mL

Detection Limit

15.6 pg/mL