



User's Manual

Growth Hormone (rat) EIA Kit



DEIA5049



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

CD's Growth Hormone (rat) EIA Kit is a competitive assay that can be used for measurement of rat GH in cell culture media or in plasma samples.

Storage

Store kit at 2-8°C or -20°C upon arrival up to the expiration date. Avoid multiple freeze/thaw cycles. For more detailed information, please download the following document on our website.

Precision

Intra-and Inter-assay coefficients of variation are 4% (n=24) and 14% (n=9).

Sensitivity

0.5 ng/mL

