



User's Manual

# Prostaglandin D2-MOX Express ELISA Kit



DEIA5054



96T, 480T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

This assay is based on the competition between PGD2-MOX and a PGD2-MOX- acetylcholinesterase (AChE) conjugate (PGD2-MOX Tracer) for a limited number of PGD2-MOX-specific rabbit antiserum binding sites.

### Storage

This kit will perform as specified if stored as directed at -20°C and used before the expiration date indicated on the outside of the box. For more detailed information, please download the following document on our website.