



User's Manual

His-Express Detection EIA Kit



DEIA5059



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

CD's His-Express Detection EIA Kit is a competitive assay designed for the rapid, semi-quantitative screening of cell lysates and affinity column fractions for His-tagged proteins.

It can be used as a substitute for SDS-PAGE to expedite the screening of affinity column fractions within a few hours.

Storage

This kit will perform as specified if stored at -20°C and used before the expiration date indicated on the outside of the box. For more detailed information, please download the following document on our website.