



User's Manual

Human Interleukin-1 β EIA Kit



DEIA5066



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

CD's IL-1 β (Human) EIA Kit is an immunometric (i.e., sandwich) EIA that permits IL-1 β measurements within the range of 0-250 pg/mL, typically with a limit of detection of 3.9 pg/mL.

Storage

Store the kit at 2-8°C. For more detailed information, please download the following document on our website.

Detection Range

0-250 pg/mL

Detection Limit

3.9 pg/mL

Sensitivity

3.9 pg/mL