



User's Manual

Leptin Receptor (Human) ELISA Kit



DEIA5071



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

This Enzyme Immunometric Assay (EIA) is based on a double-antibody sandwich technique.

Storage

Store the kit at 2-8°C. For more detailed information, please download the following document on our website.

Sensitivity

0.4 ng/mL