



User's Manual

Human FABP4 EIA Kit

REF

DEIA5073



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

CD's FABP4 (human) EIA Kit is an immunometric assay which can be used to measure FABP4 in human samples such as serum, plasma and tissue extract and tissue culture supernatant.

Storage

For more detailed information, please download the following document on our website.

Detection Limit

0.1 ng/mL