



User's Manual

Prostaglandin D Synthase (lipocalin-type) ELISA Kit

REF

DEIA5076



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

CD's L-PGDS EIA Kit is an immunometric (i.e., sandwich) EIA. The standard curve ranges from 0-100 ng/mL, with a limit of detection of 6 ng/mL. Inter- and intra-assay CV's of less than 15% may be achieved at most concentrations. L-PGDS is present in a variety of body fluids including cerebrospinal fluid, seminal fluid, and plasma. This assay has been validated using cerebrospinal fluid which contains approximately 12-30 µg/ml of L-PGDS.

Storage

This kit will perform as specified if stored as directed and used before the expiration date indicated on the outside of the box. For more detailed information, please download the following document on our website.