



User's Manual

Human Oestradiol ELISA Kit



DEIA525



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

Human Oestradiol ELISA Kit is an Enzyme Immunoassay (EIA) for the quantitative determination of total Oestradiol (E2) in human serum or plasma. For professional use only.

Storage

Store the kit reagents at 2-8°C. Do not expose to direct sunlight. DO NOT FREEZE ANY OF THE REAGENTS as this will cause irreversible damage. For more detailed information, please download the following document on our website.



Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)



Fax: 1-631-938-8221



Email: info@creative-diagnostics.com