



## User's Manual

# Human T4 ELISA Kit



DEIA528



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

Human T4 ELISA Kit is an Enzyme Immunoassay (EIA) for the quantitative determination of Thyroxine in human serum. For professional use only.

### Storage

Store components of the kit at 2-8°C upon arrival up to the expiration date. For more detailed information, please download the following document on our website.