



User's Manual

Human Free T3 ELISA Kit

REF DEIA531

 96T



RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

Human Free T3 ELISA Kit is an Enzyme Immunoassay (EIA) for the quantitative determination of Free Triiodothyronine (fT3) in human serum. For professional use only.

Storage

Store components of the kit at 2-8°C upon arrival up to the expiration date. For more detailed information, please download the following document on our website.