



User's Manual

Chlamydia Pneumoniae IgA ELISA Kit



DEIA557



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

CD' Chlamydia pneumoniae IgA test is developed for the detection of IgA antibodies specific to Chlamydia pneumoniae in human serum or plasma.

Storage

The kit should be stored at 2-8°C. For more detailed information, please download the following document on our website.