



## User's Manual

# Bovine LALBA ELISA Set



DEIA606



1000T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

Enzyme linked immunosorbent assay (ELISA) for the detection of Bovine Alpha-Lactalbumin in milk.

### Storage

Store kit at 2-8°C or -20°C upon arrival up to the expiration date. Avoid multiple freeze/thaw cycles. For more detailed information, please download the following document on our website.