



User's Manual

DNase Activity IgG ELISA Kit



DEIA6105



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

DNase Activity is a solid phase enzyme immunoassay (ELISA) for the quantitative screening of DNase Activity in human serum or EDTA-plasma.

Storage

1. Store the kit at 2-8°C.
2. Keep microplate wells sealed in a dry bag with desiccants.
3. The reagents are stable until expiration of the kit.
4. Do not expose test reagents to heat, sun or strong light during storage and usage.
5. Diluted wash buffer is stable for at least 30 days when stored at 2-8°C. For more detailed information, please download the following document on our website.

Sensitivity

0.2 mKuU/mL

