



User's Manual

M.hyop Ab ELISA Kit

REF

DEIA6123



4 x 96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The kit uses a single well blocking immunoenzymatic technique for the detection of antibodies to a membrane protein of *Mycoplasma hyopneumoniae*.

Storage

Note: Kit and diluted reagents should be stored at 2-8°C.

Samples should be stored as follows: Plasma or serum: Cold(+5°C)-max. 7days;Freeze(-20°C) -Yes;Lab Temperature(+20°C)-No For more detailed information, please download the following document on our website.