



User's Manual

Chikungunya Virus IgG capture ELISA Kit



DEIA6128



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The Chikungunya IgG-ELISA is intended for the qualitative determination of IgG class antibodies to Chikungunya in human serum.

Storage

The reagents are stable up to the expiry date stated on the label when stored at 2-8°C. For more detailed information, please download the following document on our website.

Sensitivity

The diagnostic sensitivity is >90%.