



## User's Manual

# Borrelia IgG ELISA Kit



DEIA6129



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

The recombinant *Borrelia burgdorferi* IgM-ELISA is intended for the qualitative determination of IgM class antibodies against *Borrelia burgdorferi* in human serum or plasma (citrate) and CSF (cerebrospinal fluid) offering increased diagnostic specificity and sensitivity by employing immunodominant antigens.

### Storage

The reagents are stable up to the expiry date stated on the label when stored at 2-8°C. For more detailed information, please download the following document on our website.

### Sensitivity

93%