



User's Manual

Free T4 ELISA Kit



DEIA6132



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

Competitive immunoenzymatic colorimetric method for quantitative determination of free T4 (Thyroxine) in serum or plasma.

Storage

The reagents are stable up to the expiry date stated on the label when stored at 2-8°C in the dark. Keep microtiter plate in a sealed bag with desiccant to minimize exposure to damp air. For more detailed information, please download the following document on our website.

Sensitivity

0.5 ng/l



Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)



Fax: 1-631-938-8221



Email: info@creative-diagnostics.com