



**User's Manual**

# Legionella pneumophila IgG/IgM ELISA Kit



**DEIA6148**



**96T**



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

**Creative Diagnostics**

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)**  **Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)**

---

## PRODUCT INFORMATION

### Intended Use

The Legionella pneumophila IgG/IgM Enzyme-Linked Immunosorbent Assay (ELISA) is intended for the qualitative detection of total antibodies (IgG and IgM) to Legionella pneumophila (Legionella) serogroups 1-6 in serum from patients with clinical suspicion of Legionella Disease.

### Storage

Store unopened kit at 2-8°C. For more detailed information, please download the following document on our website.



**Tel:** 1-631-624-4882 (USA)

**Tel:** 44-161-818-6441 (Europe)



**Fax:** 1-631-938-8221



**Email:** [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)