



User's Manual

Mouse H-FABP ELISA Kit (High Sensitivity)



DEIA6149



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

High sensitivity mouse H-FABP kit is offered as a tool for investigation of heart damage in mouse models of cardiovascular disease.

Storage

The reference standard stock provided with the kit should be stored at or below -20°C on receipt. The remainder of the kit should be stored at 2-8°C and the microtiter plate should be kept in a sealed bag with desiccant to minimize exposure to damp air. Test kits will remain stable until the expiration date shown on the kit package, provided that the components are stored as described above.