



User's Manual

Rat Skeletal Muscle Troponin-I ELISA Kit



DEIA6153



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

Rat Skeletal Muscle Troponin-I (SKM-TNI) ELISA Kit is for the determination of rat Troponin-I.

Storage

Store the SkM-TnI Stock vials at or below -20°C

Store the remainder of the kit at 2-8°C

Keep the microtiter plate in a sealed bag with desiccant to minimize exposure to damp air



Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)



Fax: 1-631-938-8221



Email: info@creative-diagnostics.com