



User's Manual

High Sensitivity Mouse Cardiac Troponin-I ELISA Kit

REF

DEIA6156



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

High Sensitivity Mouse Cardiac Troponin-I ELISA Kit is High Sensitivity ELISA for determination of Cardiac Troponin-I in mouse serum.

Storage

On receipt, store the lyophilized standard at or below minus 20°C. Store the remainder of the kit at 2-8°C. Keep the microtiter plate in a sealed bag with desiccant to minimize exposure to damp air.