



User's Manual

## Amyloid $\beta$ (1-38) (FL) ELISA Kit



DEIA6176



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

This IBL's assay kit is capable for the quantitative determination human amyloid  $\beta$  (1-38) which held N terminal side completely in cerebrospinal fluids and cell culture supernatant.

### Storage

Store kit reagents at 2-8°C. For more detailed information, please download the following document on our website.

### Detection Range

9.38-600 pg/mL  
(2.27-145.23 pmol/L, as molecular weight of A $\beta$  (1-38) is 4131.5)

### Sensitivity

1.28 pg/mL