



User's Manual

Amyloid- β (1-42) (FL) ELISA Kit



DEIA6178



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The kit is capable for the quantitative determination Human A β (1-42) in EDTA plasma cerebrospinal fluid and cell culture supernatant.

Storage

Store kit reagents at 2-8°C. For more detailed information, please download the following document on our website.

Detection Range

1.56-100 pg/mL(0.35-22.17 pmol/mL, as molecular weight of A β (1-42) is 4,510)

Sensitivity

0.29 pg/mL

