



User's Manual

TSH ELISA Kit



DEIA6180



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The TSH kit is intended for the quantitative determination of thyroid-stimulating hormone (TSH) in dry spots of blood of newborns.

Storage

Store all reagents at 2-8°C. For more detailed information, please download the following document on our website.

Sensitivity

3 µIU/mL



Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)



Fax: 1-631-938-8221



Email: info@creative-diagnostics.com