



User's Manual

Amyloid-beta Oligomers ELISA Kit



DEIA6190



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

This ELISA kit can detect human A β molecules which bind to anti-human A β (N) (82E1) mouse MoAb with two or more epitopes (oligomers, polymers with a protein etc.) in serum, EDTA-plasma and brain tissue extract.

Storage

Storage Condition: 2-8°C. For more detailed information, please download the following document on our website.

Detection Range

18.98-1,215 pmol/L

Sensitivity

4.41 pmol/L