



User's Manual

Resistin ELISA Kit



DEIA6200



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Human Resistin ELISA is a sandwich enzyme immunoassay for the quantitative measurement of human resistin.

Storage

Store the complete kit at 2-8°C. Under these conditions, the kit is stable until the expiration date (see label on the box). For more detailed information, please download the following document on our website.

Detection Range

1-50 ng/mL

Sensitivity

0.033 ng/mL