



User's Manual

SSA (Serum Amyloid A) ELISA Kit



DEIA6203



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The Serum Amyloid A ELISA is to be used for the quantitative determination of Human SAA in serum, plasma, buffered solution and tissue culture medium. The assay will recognize both natural and recombinant Human SAA.

Storage

Store all reagents at 2-8°C. For more detailed information, please download the following document on our website.