



User's Manual

Mouse/Rat Amyloid- β (1-42) ELISA Kit



DEIA6224



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The kit is capable for the quantitative determination Mouse/Rat A β (1-42) in serum or EDTA plasma.

Storage

Storage Condition: 2-8°C; The term of validity: 12 months

Precision

Intra-assay: 3.4%-5.8%

Inter-assay: 2.8%-8.0%

Detection Range

1.56-100 pg/mL

0.35-22.6 pmol/L (as molecular weight of A β (1-42) is 4418.0)

Sensitivity

0.05 pg/mL