



User's Manual

S 100 ELISA Kit

REF

DEIA6228



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The kit is intended for the quantitative determination of S100B (S100A1B + S100BB) in serum.

Storage

Store the kit at 2-8°C. Do not freeze. For more detailed information, please download the following document on our website.

Detection Range

The kit measures concentrations between 10 and 3500 ng/L. If S100B concentrations above the measuring range are to be expected, it is recommended to dilute samples with normal human serum prior to analysis. NOTE: The serum used for dilution should also be measured in order to determine the endogenous S100B concentration.