



User's Manual

13,14-dihydro-15-keto Prostaglandin F2 α ELISA Kit



DEIA6234



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The 13,14-dihydro-15-keto Prostaglandin F2 α EIA kit is designed for the quantitative determination of 13,14-dihydro-15-keto Prostaglandin F2 α in plasma and other sample matrices.

Storage

This kit will perform as specified if stored as directed at -20°C and used before the expiration date indicated on the outside of the box. For more detailed information, please download the following document on our website.

Detection Range

2.3-5000 pg/mL