



## User's Manual

# Rat Calcitonin Gene Related Peptide (CGRP) EIA kit



DEIA6237



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

The Rat Calcitonin Gene Related Peptide (CGRP) EIA kit is designed for the quantitative determination of rat Calcitonin Gene Related Peptide (CGRP) in serum, plasma and CSF.

### Storage

Store kit at 2-8°C or -20°C upon arrival up to the expiration date. Avoid multiple freeze/thaw cycles. For more detailed information, please download the following document on our website.