



User's Manual

Syndecan-4 ELISA Kit

REF

DEIA6243



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

This assay kit is capable for the quantitative determination Human Syndecan-4 in serum, EDTA plasma and cell culture supernatant.

Storage

Storage Condition: 2-8°C. For more detailed information, please download the following document on our website.

Detection Range

20-1,280 pg/mL

Sensitivity

3.94 pg/mL