



## User's Manual

# CrossLaps™-Serum ELISA Kit



DEIA6249



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

The Serum ELISA is an enzyme immunological test for the quantification of degradation products of C-terminal telopeptides of type I collagen in human serum and plasma.

### Storage

Store all reagents at 2-8°C and use before expiry date. When stored at 2-8°C unopened reagents will retain reactivity until expiration date. Do not use reagents beyond this date.

Unused microtiter strips must always be stored at 2-8°C in the resealable bag provided. Allow all reagents and required number of strips to reach room temperature prior to use.

Diluted, ready to use wash buffer is stable for 4 weeks when stored at room temperature (18-24°C). For more detailed information, please download the following document on our website.

### Sensitivity

0.08 nmol/L