



User's Manual

Spermatozoa-Ab/Spermatozoen-Ak ELISA Kit



DEIA6250



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

Spermatozoa-Ab/Spermatozoen-Ak ELISA is for the quantitative determination of autoantibodies against surface antigens of spermatozoa in human serum.

Storage

The kit is shipped at ambient temperature and should be stored at 2-8°C. Keep away from heat or direct sun light. The storage and stability of specimen and prepared reagents is stated in the corresponding chapters. The microtiter strips are stable up to 6 w in the broken, but tightly closed bag when stored at 2-8°C. For more detailed information, please download the following document on our website.

Sensitivity

0.14 mU/100 µL