



User's Manual

sAPP-beta wild type high sensitive ELISA Kit



DEIA6295



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

This assay kit is capable for the quantitative determination human sAPP β -w in EDTA plasma, cerebrospinal fluids and cell culture supernatant.

Storage

Storage Condition: 2-8°C. For more detailed information, please download the following document on our website.