



**User's Manual**

# **Coxiella burnetii (Q-Fever) Phase 1 IgG ELISA Kit**



**DEIA6296**



**96T**



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

**Creative Diagnostics**

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)**  **Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)**

---

## PRODUCT INFORMATION

### Intended Use

The ELISA is intended for the qualitative determination of IgG class antibodies against *Coxiella burnetii* phase 1 in the early stages of infection in human serum.

### Storage

The reagents are stable up to the expiry date stated on the label when stored at 2-8°C. For more detailed information, please download the following document on our website.



**Tel:** 1-631-624-4882 (USA)

**Tel:** 44-161-818-6441 (Europe)



**Fax:** 1-631-938-8221



**Email:** [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)