



User's Manual

Cu-Zn-Superoxid-Dismutase (SOD) ELISA Kit



DEIA6316



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The human Cu/ZnSOD ELISA is an enzyme-linked immunosorbent assay for the quantitative detection of human Cu/ZnSOD.

Storage

Store kit reagents at 2-8°C. Immediately after use remaining reagents should be returned to cold storage (2-8°C). Expiry of the kit and reagents is stated on labels.

The expiry of the kit components can only be guaranteed if the components are stored properly, and if, in case of repeated use of one component, the reagent is not contaminated by the first handling. For more detailed information, please download the following document on our website.

Detection Range

11.0-16.0 ng

Sensitivity

0.04 ng/mL