



User's Manual

[pSer473/474]Akt1/2 ELISA Kit



DEIA6323



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The [pSer473]Akt1/2 REIA kit is a complete kit for the quantitative determination of Akt1 and 2 from human, mouse, and rat cell lysates.

Storage

Upon receipt store all components at 4°C. For more detailed information, please download the following document on our website.

Detection Range

17.5-560 pg/mL

Sensitivity

5.55 pg/mL

