



User's Manual

Bradykinin ELISA Kit



DEIA6335



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Bradykinin Enzyme-Linked Immunosorbent Assay (ELISA) kit is a complete kit for the quantitative determination of Bradykinin in plasma, serum and urine.

Storage

All components of this kit, except the Standard, should be stored at 4°C upon receipt. The Standard should be stored at -20°C. Shipping conditions may not reflect storage conditions. For more detailed information, please download the following document on our website.

Sensitivity

24.8 pg/mL



Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)



Fax: 1-631-938-8221



Email: info@creative-diagnostics.com