



User's Manual

β-Catenin ELISA Kit



DEIA6338



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The β-catenin ELISA kit is a complete kit for the quantitative determination of β-catenin in cell lysates.

Storage

All components of this kit, except the β-catenin Standard, are stable at 4°C until the kit's expiration date. The β-catenin Standard must be stored at or below -20°C. For more detailed information, please download the following document on our website.

Sensitivity

33.3 pg/mL



Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)



Fax: 1-631-938-8221



Email: info@creative-diagnostics.com