



User's Manual

# Human HO-1 ELISA Kit



DEIA6384



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

The HO-1 (human), ELISA kit provides a method to detect and quantitate Heme Oxygenase (HO-1) in samples from human origins.

### Storage

All reagents are stable as supplied at 4°C, except the Recombinant HO-1 Standard, which should be stored at -20°C. For optimum storage, the Recombinant HO-1 Standard should be aliquotted into smaller portions and stored at -20°C.

Unused wells of the Anti-HO-1 Immunoassay Plate should be resealed in the foil pouch provided and stored at 4°C until the kits expiry date.

### Sensitivity

0.78-25 ng/mL