



User's Manual

13(S)-HODE ELISA Kit



DEIA6386



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The 13-(S)-HODE ELISA kit is a competitive immunoassay for the quantitative determination of free 13(S)-HODE in biological fluids.

Storage

All components of this kit, except the conjugate, are stable at 4°C until the kit's expiration date. The conjugate must be stored frozen at -20°C. For more detailed information, please download the following document on our website.

Sensitivity

1.6 ng/mL

