



## User's Manual

# HSP70B ELISA Kit



DEIA6397



96T



This product is for research use only and is not intended for diagnostic use.



For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

The HSP70B' ELISA kit provides a sensitive and specific method to detect and quantitate HSP70B' in cell lysates, tissue extracts and serum samples from human origin.

### Storage

All reagents are stable as supplied at 4°C, except the Recombinant HSP70B' Standard, which should be stored at -20°C. For optimum storage, the Recombinant HSP70B' Standard should be aliquotted into smaller portions and stored at -20°C. Avoid repeated freeze thaw cycles.

Unused wells of the Anti-HSP70B' Immunoassay Plate should be resealed with desiccant in the foil pouch provided and stored at 4°C until the kits expiry date. For more detailed information, please download the following document on our website.

### Sensitivity

0.062 ng/mL