



User's Manual

24(S)-Hydroxycholesterol ELISA Kit

REF

DEIA6398



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The 24(S)-Hydroxycholesterol (24-OHC) Enzyme-Linked Immunosorbent Assay (ELISA) kit is a complete kit for the quantitative determination of 24-OHC in tissue culture media, cerebral spinal fluid and tissue homogenate samples.

Storage

All kit components should be stored at 4°C upon receipt. Shipping conditions may not reflect storage conditions. For more detailed information, please download the following document on our website.

Sensitivity

0.78 ng/mL